



ABRASP Meeting

Sao Paulo, October 20th, 2011





Ultra-rapid Microbiology using ChemScan[®] RDI & Regulatory Requirements for Validation of the RMM



Solid Phase Cytometry ChemScan[®] RDI



Content:

1. Principle of the Laser Scanning Cytometry Technology
2. Applications
3. An example of implementation of ChemScan[®] RDI




Principle of the Laser Scanning Cytometry Technology


ChemScan[®] RDI



Chemunex technology avoiding the need for cell growth



Fluorassure[®] Viability Markers



- 1 Accumulation of the viability substrate in the cell
Membrane integrity

- 2 Activation of the viability substrate by the enzyme in the cytoplasm
Enzyme activity


ChemScan analysis : A simple three step procedure




1. Filtration



2. Cell Labelling




3. Laser scanning



1. Sample Filtration



Large volumes can be tested using :
standard filtration units or ready to use FIFUs



0.4µm polyester track-etched membranes (= ChemFilter)

2. Analysis Protocol with FIFU



1. FILTRATION

2. LABELLING


3. LASER SCANNING





Sample

Filtration


Activation



Labelling



Chem Scan Analysis




Pad saturated with activation solution



Pad saturated with labelling solution

3. Laser Scanning






→ 3 ou 5 min scanning for the all membrane

Fluorescent cell labelling




Incident light

488 nm




Fluorescent Signal




Fluorescence Based Cell Labels

3. Scanning is fully overlapping




3. Data Processing






Signal Shape



Labelled Micro-organism




Autofluorescent Particle




Signal Discrimination

Data Map = Total count



Results Map = labelled microorganisms



Rejected Background

- Autofluorescent Particles
- Membrane Fluorescence
- Electronic Noise

Results Display



Chemsan RDI Session = 31H09 - [Log Window]

File Microscope Users Window Help

Session ID: 31H09 Operator ID: USER Session Date: 31/08/09 Session Time: 13:49:57 Session Comment:

Scan Name	Scan Date	Scan Time	Sample ID	Batch ID	Comment	Results	Microscope Results	Scan User ID	Scan Status
31H09000	31/08/09	13:51:59	Standard C3	C903307		711	86	USER	SCAN VALIDATED
31H09001	31/08/09	14:01:14	water sample	001		8		USER	READY TO VALIDATE

Chemsan RDI Session = 31H09 - [Scan Maps - d:\live\31H09\31H09001]

File View Users Window

Cursor:
 X (mm): 00.000
 Y (mm): 12.618
 Results
 Nearest Spot: 8

Mouse Click Options:
 Zoom Mode:
 Nearest Spot Mode:
 Next Spot Mode:
 Absolute Mode:
 Move Microscope:

Begin Validation
 Microscope Setup
 Print
 Close


Validation:

Bacterium	+	6	Load
Unobserved	⊖	0	Save
Particle	■	2	Quit
Total:		8	


Results

Current Spot: 8
 Spots: 8


Microscope Validation



1. Membrane holder



2. Automated Microscope Stage




3. Validation of the Scan map


Microscope objective

Monitized platform


Labelled Microorganisms



Candida



Bacillus



Mould

Direct viability labelling demonstrated with wide range of microorganisms

Bacteria Gram -

Achromobacter xylosoxydans
Aeromonas hydrophila
Agrobacterium radiobacter
Alcaligenes eutrophus
Alcaligenes faecalis
Burkholderia cepacia
Burkholderia diminuta
Burkholderia pickettii
Caulobacter sp.
Cedecea lapagei
Citrobacter diversus
Citrobacter freundii
Comamonas terrigena
Edwardsiella hoshinae
Enterobacter aerogenes
Enterobacter agglomerans
Enterobacter cloacae
Enterobacter gergoviae
Enterobacter sakasaki
Enterobacter intermedium
Erwinia Sp.
Escherichia coli
Escherichia coli HB 101
Escherichia coli 0126 :B16
Flavobacterium Sp.
Klebsiella oxytoca
Klebsiella planticola
Klebsiella planticola
Klebsiella pneumoniae
Klebsiella terrigena

Kluyvera Sp.
Moraxella sp.
Pasteurella aerogenes
Proteus mirabilis
Pseudomonas diminuta
Pseudomonas aeruginosa
Pseudomonas alkagenèse
Pseudomonas mesophilica
Pseudomonas putida
Pseudomonas fluorescens
Pseudomonas stutzeri
Salmonella choleraesuis
Salmonella indiana
Salmonella typhimurium
Salmonella eboni
Salmonella sp.
Salmonella virchow
Serratia marcescens
Shigella sonnei
Xanthomonas maltophilia
Yersinia enterocolitica

Bacteria Gram +

Aerococcus viridans
Bacillus anthracis
Bacillus amyloliquefaciens
Bacillus cereus
Bacillus circulans
Bacillus coagulans
Bacillus globigii

Bacillus lentus
Bacillus licheniformis
Bacillus megaterium
Bacillus mycoides
Bacillus pumilus
Bacillus sphaericus
Bacillus stearothermophilus
Bacillus subtilis
Bacillus thuringiensis
Bacteroides fragilis
Bacteroides thetaiotamicron
Bacteroides vulgatus
Clostridium acetobutylicum
Clostridium bifermentans
Clostridium butyricum
Clostridium perfringens
Clostridium sporogenes
Clostridium tyrobutyricum
Corynebacterium aquaticum
Corynebacterium pseudodiphtheriticum
Enterococcus faecium
Enterococcus faecalis
Fusobacterium nucleatum
Lactobacillus acidophilus
Lactobacillus brevis
Lactobacillus buchneri
Lactobacillus bulgaricus
Lactobacillus casei casei
Lactobacillus casei
Lactobacillus cellobiosus
Lactobacillus curvatus
Lactobacillus delbrueckii

Lactobacillus fermentum
Lactobacillus leichmannii
Lactobacillus plantarum
Lactobacillus lactis
Lactobacillus sake
Lactobacillus sp.
Leuconostoc oenos
Leuconostoc Sp.
Listeria innocua
Listeria monocytogenes
Micrococcus luteus
Mycobacterium bovis
Mycobacterium parafortuitum
Mycobacterium megmatris
Mycobacterium tuberculosis oerskovia sp.
Pediococcus damnosus
Pediococcus pentosaceus
Porphyromonas canoris
Porphyromonas gingivalis
Propionibacterium acnes
Staphylococcus aureus
Staphylococcus epidermidis
Staphylococcus hominis
Staphylococcus warneri
Staphylococcus xylosus
Streptococcus faecalis
Streptococcus salivarius
Streptococcus thermophilus
Streptococcus viridans
Thiobacillus ferrooxidans

Direct viability labelling demonstrated with wide range of microorganisms



Yeast

Acremonium kiliense
Candida albicans
Candida ciferii
Candida colliculosa
Candida famata
Candida famata
Candida fumentans
Candida humicola
Candida humicola
Candida krusei
Candida luxitaniae
Candida magnolia
Candida parapsilosis
Candida pelliculosa
Candida tropicalis
Cryptococcus albidus
Debaryomyces hansenii
Galactomyces geotrichum
Geotrichum candidum
Hansenulasporea uvarum
Hansenula anomala
Kloeckera japonica
Kloeckera Apis apiculata
Pichia anomala
Pichia guilliermondii
Pichia menbrena faciens
Rhodotorula rubra
Saccharomyces baillii
Saccharomyces bisparus
Saccharomyces cerevisiae
Saccharomyces rosei

Torulopsis candida *Torulopsis*
inconspicua *Torulopsis maris*
Torulosporea delbrueckii
Zygosaccharomyces baillii
Zygosaccharomyces rouxii

Mould

Acremonium Sp.
Aspergillus versicolor
Aspergillus versicolor
Aspergillus fumigatus
Aspergillus niger
Basydiomycetes Sp.
Bassochlamys fulva
Byssochlamys Sp.
Cladosporium cladosporioides
Epicocum nigrum ou altenaria
Fusarium oxysporum
Fusarium oxysporum
Fusarium gramineatum roseum
Humicola fuscoatra
Mucor circinelloides
Mucor plumbeus
Mucor racemosus
Mucor Sp.
Neosartoea Sp.
Penicillium decumbens
Penicillium expansum
Penicillium frequentans
Penicillium roquefortii
Rhizopus Sp.
Rhodotorula rubra
Rhizopus oligosporus
Scopulariopsis candida
Trichoderma Sp.



ChemScan[®] RDI Applications



Pharmaceutical applications



Exemple of applications:

Sterility Test

Sterility Test for filterable products



Less than 4 hours

Biourden

TVC Bioburden for in-process

TVC Bioburden for raw material

TVC Bioburden for end-product

90 min



90 min

< 3 hours

Envirommental controls

TVC Bioburden for pharmaceutical water

Air monitoring using Coriolis

Surface monitoring using ChemSwab

90 min



< 3 hours

< 3 hours



Biotechnology

Contaminations of cell cultures

< 2 hours

Control of fermentations

90 min

1. Pharmaceutical water testing



- Water Testing demands a high sensitivity & enumeration of microorganisms (quantitative results)
- Current growth based methods have limitations :
 - can lead to underestimation since certain micro-organisms, such as spores, stressed cells & fastidious cells are unable to grow on TSA agar
 - Incubation period between 48-72 hours leads to delay in detection of problems and intervention
 - retrospective to the use of water in the production process (2-14 days to result)
 - can be a compromise between speed & sensitivity
 - delay investigative & corrective activities

1. Pharmaceutical water testing



❖ Applications evaluated :

- Characterisation of new and existing water systems
- Biofilm monitoring
- Routine water system trending

❖ Acknowledgements :

- Gunter Gapp - Novartis
- Sylvie Guyomard – Sanofi Aventis
- Pascale Nabet – Sanofi Aventis
- Jean Scouart - UCB Pharma


Ref : Evaluation of the applications of a system for real time microbial analysis of pharmaceutical water systems, EJPS; 1999, 4(4): 131-136

1. Pharmaceutical water testing



Case study :


- two purified water systems analysed :
 - 1st: newly constructed system designed to minimise opportunities for biofilm formation
 - 2nd: older system prone to biofilm formation and required regular attention
- both systems analysed with the ChemScan[®] RDI and using R2A plates incubated for 5 days at 32°C



1. Pharmaceutical water testing



Monitoring of 'new' water system




- all results in the range of 0-10 counts / ml
- good correlation between ChemScan & plate counts

1. Pharmaceutical water testing



Monitoring of 'old' water system



- counts higher than seen in the 'new' water system by both methods
- higher counts on ChemScan compared to plate count data
- related to the ability of ChemScan to detect stressed cells, spores and fastidious microorganisms not recovered by high nutrient plate count methods (e.g. TSA)
- plate count depends on the ability of the organism present to grow under the particular culture condition used

1. Pharmaceutical water testing



Conclusions:

- ➔ ChemScan typically **equivalent** to plate count methods in water systems which are in control
 - results within current action levels

- ➔ System has substantially **higher sensitivity** for detection of stressed cells, spores and fastidious organisms
 - sensitivity in 90 minutes at least as good as with low nutrient media incubated for extended periods (5-14 days)
 - improved detection of biofilm formation
 - opportunities for early detection of contamination not seen with current plate count methods (clear financial benefits)
 - opportunity for routine real time monitoring of water systems without compromising sensitivity

2. Bioburden in Intermediates




Moving microbial controls from Product to Process :

- **Microbial quality of incoming raw materials**
- **Bioburden on bulk solutions**

Fast analysis of bulk products before filling

Faster turn-over of production tanks

Fast detection, localisation & correction of contamination



2. Bioburden in Intermediates




Pharmaceutical Products (examples):

- ✓ Analgesics
- ✓ Antibiotics
- ✓ Contact lens washing solutions
- ✓ Detergents and cleaning products
- ✓ Antiseptic solution
- ✓ Heparin
- ✓ Nasal solutions
- ✓ Peritoneal dialysis solution
- ✓ Shampoo
- ✓ Sugar solutions
- ✓ Vaccines
- ✓ Vitamins

3. Surface Monitoring using ChemSwab[®]





Quantitative rapid surface monitoring using ChemSwabs (Flocked Fiber swabs) and ChemScan[®] RDI rapid method



+



3. Surface Monitoring using ChemSwab®



SWAB RECOVERY CAPACITY

Recovery from flocked swabs up to **3 TIMES HIGHER** than traditional swabs (**60% versus 20%**).

SWAB RELEASE CAPACITY

Flocked swabs release **4 TIMES HIGHER** than traditional swabs, (**90% versus 20%**).

Ref: GSK Article in PDA Journal

3. Surface Monitoring using ChemSwab®



Comparative Study

Material and Methods :

- 50uL of the suspension on the swabs
- Filtrate solution on ChemFilter CB04 (with a funnel) also compatible with FIFU
- Rinse the Cone with a specific solution / filtrate the rinsing volume
- Analyze the sample with "Surface monitoring using ChemSwabs" protocol on ChemScan® RDI

Tested strains :

Staphylococcus epidermidis SE
Bacillus subtilis Spores BSS
Escherichia coli EC
Aspergillus brasiliensis spores ABS
Candida albicans CA
Burkholderia cepacia BC

ChemSwab®



3. Surface Monitoring using ChemSwab[®]



Comparative Results

(CFU = control plate, Scan = ChemSwabs + Scan RDI protocol
(number of micro organism detected)) :

Strains	EC		CA		ABS		BSS		SE		BC	
	Scan	CFU	Scan	CFU	Scan	CFU	Scan	CFU	Scan	CFU	Scan	CFU
	91	81	158	103	79	66	141	91	158	156	105	103
	95	83	176	98	92	61	148	77	176	159	89	78
	81	87	143	99	91	73	172	86	143	155	68	92
	-	-	170	102	136	80	125	99	170	168	84	93
Mean	91	84	162	101	100	70	147	87	162	160	87	92
Correlation Scan/CFU	108%		160%		143%		169%		101%		95%	

4. Air Monitoring using Coriolis[®]μ



Portable Air Sampler for bio-aerosol collection

Your results in few hours combining



Coriolis[®]μ

with ChemScan[®]RDI



Designed and developed by:




Designed and developed by:

4. Air Monitoring using Coriolis[®]_μ



Operating principle of the technology



Before collection

Cone ready to use
with the adapted
collection liquid

During collection

1. Vortex beginning
Air and particles
come into the cone
forming the vortex

2. During the vortex
Aspirated particles
are centrifuged
with the liquid on the
cone's inner surface

After collection

Particles
concentrated
in the liquid





4. Air Monitoring using Coriolis[®] μ



Comparative Study

Material and Methods :

- 1,5m³ sampling: 5 minutes at 300 L/min with Coriolis[®] μ using collection liquid
- Filtrate solution on ChemFilter CB04 (tests with a funnel) also compatible with FIFU
- Rinse the Cone with a specific solution / filtrate the rinsing volume
- Analyze the sample with "Air monitoring using Coriolis" protocol on ChemScan RDI

Tested strains :

Staphylococcus epidermidis SE
Bacillus subtilis Spores BSS
Escherichia coli EC
Aspergillus brasiliensis spores ABS
Candida albicans CA
Burkholderia cepacia BC
Micrococcus luteus ML



4. Air Monitoring using Coriolis[®]_μ



Comparative results

(CFU = control plate, Scan = Coriolis + Scan RDI protocol
(number of microorganism detected)) :

Strains	EC		CA		ABS		ML		BSS		SE		BC	
	Scan	CFU	Scan	CFU	Scan	CFU	Scan	CFU	Scan	CFU	Scan	CFU	Scan	CFU
	55	52	91	97	57	45	101	89	72	98	91	33	64	81
	64	56	69	78	45	46	113	85	69	91	74	48	59	67
	68	55	94	82	72	57	88	77	93	88	80	41	62	69
Mean	62	54	85	86	58	49	101	84	78	92	82	41	62	72
Correlation Scan/CFU	115%		99%		118%		120%		85%		201%		85%	

Evaluation of Coriolis® microbial air sampler coupled with RMM

Alternative solution for rapid airborne contamination control



Abstract

Environmental contamination control in cleanrooms, Bertin Technologies (France) has developed a technology dedicated to the monitoring of airborne particles. The goal is to propose a **sampling method compatible with Rapid Microbiological Methods (RMM)** in order to get reliable & specific results in a few minutes. The impact on pharmaceutical production of time-to results within impaction method, the Coriolis® μ air sampler has been validated according to ISO14698-1 in terms of physical and biological efficiencies (HPA study – July 2008) : the Coriolis® μ is as efficient as the traditional method and even better for high particles diameter (results available on www.coriolis-technologies.com). Pharmunex has also validate a protocol coupling Coriolis® μ with ScanRDI® (cytometry) and allows to get the results in only 3 hours (from sampling to results). This study aims at completing this data and at testing different RMMs on the samples collected with Coriolis® μ.

Context



Microbiological quality control of environments aims at ensuring the **quality of products** in case of cleanrooms production / **health and safety** of workers and exposed people. Reliable measurements of microbial contamination depend on the choice of an adapted air sampler / a representative sample from the environment / the limitation of losses due to a failure of the sampler to capture particles containing micro-organisms or to due to the formation of viable micro-organisms during collection so that formation of visible colonies on agar surfaces will not occur. **Objective** : realize the feasibility study of the Coriolis® μ air sampler coupled with RMM (Scan RDI) in order to implement it as a **new solution for rapid investigation** in case of contamination and for **monitoring in production sites**.


Material & methods

Sampling has been carried out in one of the Micro Lab using either traditional air sampler (impactor) or Coriolis® μ. Samples collected with the Coriolis® μ were further processed by filter-plate or ScanRDI method (RMM).

COMPARISONS = Impactor vs Coriolis® μ (Bertin Technologies)
 Coriolis® μ air sampler is based on a patented **cyclonic technology**: it concentrates airborne particles into liquid media.



SCANRDI METHOD = ScanRDI (AES CHEMUNEX)
 ScanRDI is based on cytometry and uses a double discrimination key: viability and cell membrane integrity.



Keywords : biocontamination – airborne particles – microbial air

Results

- Experiment A** = comparison of Coriolis® μ and impactor (table 1)
 - Different areas controlled, equipments placed side by side, 2 replicates
 - Filtration of the whole sample from Coriolis® sampling (0,45μm filter) + filter placed on TSA agar plate
 - Incubation of plates 20-25°C for 72-96h + 30-35°C for 48-72h
- Experiment B** = Coriolis® μ sampling and ScanRDI detection (table 2)
 - Each liquid sample from Coriolis® splitted in A&B aliquots
 A = filter plate method
 B = filtered and analysed with ScanRDI

Table 1 - the impactor and Coriolis® μ give equivalent results on the air samples in the lab

Samples	Impaction	Coriolis® μ + filter plate
Biohood	0	1 M
Biohood	0	0
Counter top 2	6 (4 M + 2 B)	2 B
Counter top 2	4 (2 M + 2 B)	1 M
Counter top 3	6 (3 M + 3 B)	8 (7 M + 1 B)
Counter top 3	2 B	6 (3 M + 3 B)

M = Mold - B = Bacteria

Table 2 - ScanRDI analysis gives higher counts than the filter plate method does, especially for contaminated air samples

Samples	Coriolis® + Scan RDI	Coriolis® + filter plate
Biohood	6	0
Biohood+hands	1	0
Counter top 3	9	2 (1 M + 1 B)
Counter top + walking	66	7 (3 M + 4 B)
Counter top + hands	103	8 B

These higher results are partly due to the **viable but non-culturable (VNC)** microbials present in the air which can not be detected by the traditional impaction method although they can be pathogenic.

Discussion

Liquid sample → Access to alternative methods → Rapid results (RMM beyond cultivable flora)

The Coriolis® system is an easy to use portable air sampler that collects air samples into liquid media. This **allows quick and reliable detection when coupled with RMM technologies** as demonstrated into this study.

The feasibility study conducted in Saint-Louis Micro Lab indicates that Coriolis® system collects airborne microorganisms in a comparable amount to the impaction system does. ScanRDI can be used as Rapid Microbiology Method coupled with the Coriolis® system in order to give rapid results (around 3 hours from sampling to result). It is also shown here that the amount of microorganisms detected with the ScanRDI is higher than with the traditional method as far as it is not based on cultivability. Appropriate alert and action limits may thus need to be re-evaluated for ScanRDI results. This couple of innovative technologies fits for the Environmental Monitoring application and could be implemented for **investigation and routine monitoring** in production sites and in critical areas and **cleanroom environments**.

Acknowledgement



• Pfizer study Team: Dr Lin Chen (responsible of the study), S. Fennel
 • Bertin Technologies: Alexandra Guerin, Quiterrie Desjonquères
 • AES CHEMUNEX: Pierre Barbez



Conclusion :
an exemple of
implementation
of the ChemScan[®] RDI




How useful is microbiological testing?



How useful is microbiological testing of end product ?

Lot: 0.1% defectives → 10 samples analyzed: Probability of detection ~ 1%

Lot: 20% defectives → 10 samples analyzed: Probability of detection ~ 35%



How Can we move Forward ?



- Improving sampling sensitivity
 - Air and gases
 - Surfaces and personnel

- Implementation of Rapid Micro Technology
 - Selecting the suitable for the purpose
 - Defining an appropriate strategy

- Moving microbial controls from Product to Process
 - Microbial quality of incoming materials : Raw material testing and In-process Bioburden
 - Environmental monitoring : Air and surface monitoring
 - Water for pharmaceutical use

Recent Case of implementation for PAT



- Non sterile nasal product
- Deployment of environmental monitoring, intermediates and waters data for Real Time Release (RTR) of drug products manufactured under conventional aseptic process
- 6 critical points from risk analysis

Recent Case of implementation for PAT



From which:

➤ Process Water	→	ChemScan [®] RDI - 90 min
➤ Bulk Product	→	ChemScan [®] RDI - 3 hours
➤ Surface	→	ChemScan [®] RDI - 3 hours
➤ Air	→	Other RMM - 24 hours
➤ Filled Vials	→	ChemScan [®] RDI - 3 hours

This manufacturing site receives FDA inspection and formal approval in July 08 to use ChemScan[®] RDI system as a part of a rapid microbiological in-process monitoring of a non sterile product **eliminating the need for end product microbial testing prior to release.**



ChemScan[®] RDI Validation



Alternative methods validation



3 standards on the validation of alternative methods:

- European Pharmacopoeia, chapter 5.1.6 “Alternative Methods for Control of Microbiological Quality”
 - Description of some alternative methods in which the ChemScan[®] RDI description (solid phase cytometry) and BactiFlow[®] ALS description (flow cytometry)
- USP, chapter 1223 “Validation of Alternative microbiological methods”
- PDA, Technical Report 33, “The Evaluation, Validation and Implementation of New Microbiological Methods”

Validation steps



1. System Qualification
 - FAT (Factory Acceptance Test)

2. Qualification of the installation of the system and its main operational functions
 - Commissioning
 - IQ : Installation Qualification
 - OQ : Operational Qualification

3. Performance Qualification of the system
 - PQ1 : Performance Qualification 1

4. Performance Qualification in real conditions (with samples)
 - PQ2 : Performance Qualification 2

Performance Qualification



Validation Parameter	Qualitative Test		Quantitative Test	
	European Pharmacopeia	US Pharmacopoeia	European Pharmacopeia	US Pharmacopoeia
Accuracy	Yes	No	Yes	
Linearity	No		Yes	
Precision	Yes	No	Yes	
Limit of Detection	Yes	Yes	No	Yes
Limit of quantification	No		Yes	
Assay Range	No		Yes	
Specificity	Yes		Yes	
Robustness / Ruggedness	Yes		Yes	

Validation of the ChemScan[®] RDI



ChemScan[®] RDI

Validation of the ChemScan[®] RDI with quantitative applications



- Commissioning

1. Installation Qualification

2. Operational Qualification

3. Performance Qualification 1

4. Performance Qualification 2

Customer
service

3 days

Customer
service

1 week

Or

Or

Laboratory

Several weeks

Laboratory

5 – 6 weeks*

Depend on
the laboratory

Validation of the ChemScan[®] RDI with qualitative applications



- Commissioning

1. Installation Qualification

2. Operational Qualification

3. Performance Qualification 1

4. Performance Qualification 2

Customer
service

3 days

Customer
service

1 week

Or

Or

Laboratory

Several weeks

Laboratory

4 weeks*

Depend on
the laboratory

Doc. Ref. : 200-D0225-11	Performance Qualification 1	Page : 1 of 54
Author : L. Jost	Protocol and Report	
Date : 28 October 2008	(Chem)Scan® RDI Unit For TVC Bioburden, Fungi and Scan Bio II protocols	

PERFORMANCE QUALIFICATION 1
ANALYTICAL PERFORMANCE PROTOCOLS AND REPORT
 AES CHEMUNEX, France

(Chem)Scan®RDI UNIT
for TVC Bioburden, Fungi, Scan Bio II
protocols

Document Reference:	200-D0225-11
Date of Issue:	28 October 2008
Page:	1 of 54

Affix the self adhesive control label here

Protocol Approvals

Authors:

Name	Job Title	Company	Signature:	Date:
L. Jost	Application Specialist	AES CHEMUNEX		

Approved By:

Name	Job Title	Company	Signature:	Date:
J-L. Droocourt	Scientific & Technical Director	AES CHEMUNEX		

Report Approvals

Name	Job Title	Company	Signature:	Date:

Doc. Ref. : 200-D0225-11	Performance Qualification 1	Page : 2 of 54
Author : L. Jost	Protocol and Report	
Date : 28 October 2008	(Chem)Scan® RDI Unit For TVC Bioburden, Fungi and Scan Bio II protocols	

TABLE OF CONTENTS

1. OBJECTIVES AND SCOPE	3
2. DESCRIPTION	5
3. RESPONSIBILITIES	7
4. DATA COLLECTION AND ENTRY	8
5. MATERIALS / EQUIPMENT	9
6. OUTSTANDING ISSUES	10
7. SAFETY PRECAUTIONS AND APPROVAL TO COMMENCE WORK	11
8. IDENTIFICATION OF PERSONNEL	12
9. QUALIFICATION STRATEGY	13
10. SOP VERIFICATION	15
11. SYSTEM PERFORMANCE MONITORING	17
12. RECOMMENDED METHODS FOR STANDARD MICRO-ORGANISM SAMPLE PREPARATION	20
13. ANALYTICAL CHECKS	21
14. PREVENTATIVE MAINTENANCE	31
15. TRAINING	32
16. CHANGE CONTROL	33
17. LIST OF APPENDICES	34
18. NON-CONFORMANCE LIST	35
19. SUMMARY AND CONCLUSIONS	37
20. DOCUMENT REVIEW AND APPROVALS	38
21. FINAL APPROVALS	39
22. REFERENCE DOCUMENTS	40
23. KEY TO ABBREVIATIONS	41
24. PROTOCOL REVISION HISTORY	42

Doc. Ref. : 200-D0225-11	Performance Qualification 1	Page : 13 of 54
Author : L. Jost	Protocol and Report	
Date : 28 October 2008	(Chem)Scan® RDI Unit For TVC Bioburden, Fung and Scan Bio II protocols	

9. QUALIFICATION STRATEGY

9.1 Qualification source references

The qualification strategy as defined in the PQ1 has been compiled from the following key analytical references:

- US Pharmacopeia 31, "<1223> Validation of Alternative microbiological methods"
- PDA Technical Report 33, "The Evaluation, Validation and Implementation of New Microbiological Methods"
- European Pharmacopoeia 6th Edition, "5.1.6 Alternative Methods for Control of Microbiological Quality"

The key analytical parameters that need to be tested are for the PQ1 of a quantitative application:

- **Accuracy:** This is defined as the closeness of the test results obtained by the alternative method to the value obtained by the pharmacopoeial method. Accuracy must be demonstrated across the practical range of the test.
- **Linearity:** This is defined as the ability to produce results that are proportional to the concentration of micro-organisms present on the sample within a given range.
- **Precision:** This is defined as the degree of agreement among individual test results when the procedure is applied repeatedly to multiple samplings of homogeneous suspensions of micro-organisms under the prescribed conditions.
- **Limit of quantification :** This is defined as the lowest number of micro-organisms that can be accurately counted
- **Limit of Detection:** This is defined as the lowest number of micro-organisms in a sample that can be detected under the stated experimental conditions.
- **Assay Range:** This is defined as the interval between the upper and lower levels of micro-organisms that have been demonstrated to be determined with precision, accuracy and linearity using the method as written.
- **Specificity:** This is defined as the ability of the method to accurately detect a required range of micro-organisms that be present in the sample under test.
- **Robustness (and Ruggedness):** This is defined as the measure of its capacity to remain unaffected by small but deliberate variations in method parameters and provides an indication of its reliability under a variety of normal test conditions, such as different analysts, instruments, batches of reagents and laboratories.

Doc. Ref. : 200-D0225-11	Performance Qualification 1	Page : 22 of 54
Author : L. Jost	Protocol and Report	
Date : 28 October 2008	(Chem)Scan® RDI Unit For TVC Bioburden, Fung and Scan Bio II protocols	

Check 1- Determination of Accuracy and Linearity

OBJECTIVES

To compare the microbiological counts derived from testing samples with the (Chem)Scan® RDI with those from a traditional microbiological technique (a Reference standard) to determine Accuracy and Linearity.

The studies should include all the organisms of relevance to the production process assayed as single pure cultures. The identity of the micro-organism should be proven by a defined method and reference made as to the location of this information.

The range of analysis should be determined according to the range of interest for the routine analysis: The highest concentration of the range should be higher than the acceptance limit and the lowest concentration of the range should correspond to the limit of quantification of the technology (in most of the case, 5 cfu / filtered volume is acceptable).

The limit of detection of the technology will be calculated later on.

It is recommended to analyze five suspensions where concentrations are within the range of analysis, and repeat the assay 3 times.

N.B: Assay the dilutions in 3 series of five singletons and perform the plate assays at the same time. This removes any bias that may be caused by the cultures multiplying between assays.

Two other working sessions are then carried out under conditions of maximum variability (different reagents, different operators, different days, etc).

As the reproducibility is referring to "the use of microbiological method within the same laboratory over a short period of time using different analysts with the same equipment" (Source: PDA Technical Bulletin 33), it is recommended to change analysts for the different working sessions.

The following table shows an example for a range of interest between 5 and 100 cfu / filtered volume:

	Working day 1			Working day 2			Working day 3			
Dilution A	160 cfu	160 cfu	160 cfu	160 cfu	160 cfu	160 cfu	160 cfu	160 cfu	160 cfu	Dilution factor 2
Dilution B	80 cfu	80 cfu	80 cfu	80 cfu	80 cfu	80 cfu	80 cfu	80 cfu	80 cfu	Dilution factor 2
Dilution C	40 cfu	40 cfu	40 cfu	40 cfu	40 cfu	40 cfu	40 cfu	40 cfu	40 cfu	Dilution factor 4
Dilution D	10 cfu	10 cfu	10 cfu	10 cfu	10 cfu	10 cfu	10 cfu	10 cfu	10 cfu	Dilution factor 2
Dilution E	5 cfu	5 cfu	5 cfu	5 cfu	5 cfu	5 cfu	5 cfu	5 cfu	5 cfu	

The same data will be used to determine the Linearity, Accuracy, Precision and Limit of Quantification of the alternative method.



CHEMUNEX 
The Rapid Microbiology Company

Subsidiary of



presents

**Muito
Obrigado**

